


# REUSABLE CYSTOSCOPY

**Comfort in Every Cystoscopy**  
**Cysto-Nephroscope**



# RE-USABLE CYSTOSCOPE

## Standard equipped with a suction button

Remove small tissue samples or debris from the bladder; Remove the fluids such as urine or irrigation fluid which accumulate in the bladder, resume the operational view that is affected by the blood diffusion.

### Larger Working Channel

With a large working channel, our cystoscope delivers enhanced access and expanded options, empowering clinicians to perform procedures with ease and efficiency.

### Advanced insertion tube design

Flexible cystoscopes are less rigid and thinner than rigid cystoscopes, reducing the discomfort for patients and allowing easier access through the urethra. Additional marks are added on the insertion tube, indicating the positioning of the insertion tube.

### Insertion tube rotation

Extra rotation knob at the bottom of the cystoscope handle, providing higher maneuverability for surgeons during the procedures.

### Customizable Buttons

Two in-scope buttons that can be customized with different functions, such as white balance and image recording.

### Enhanced bending section design

180 bending angle proves sufficient  
Angle of viewing the t

### Higher Portability

Compatible with HugeMed's MS-8 image processor, with its light and compact design, it is especially convenient for out-patient department and ICU daily examination.



Model	C Y-Q5 8	C Y-Q5 2	C Y-Q4 0	C Y-Q3 0
Insertion Diameter	17.4FR	15.6FR	12.0FR	9.0 FR
Channel Diameter	9.0 FR	8.4 FR	6.0 FR	3.6 FR
Working Length	350mm	350mm	350mm	350mm
View Angel	120°	120°	120°	120°
View Depth	3-5 0mm	3-5 0mm	3-5 0mm	3-5 0mm